




WELTRAUM SOFTWARE UNINSTALLATION




Click Start->Run from windows Desktop



Type Control Panel. You May see this Page




Or




Click on Programs and Features.


Removing WELTRAUM Software



Look for WELTRAUM Software. Right Click and Click on Uninstall/Change.



Click OK to Proceed.



Click Yes to all



Once Completed Click OK

Removing SQL Server Instance.

In the Program and features Screen

Control Panel Home

Uninstall or change a program

View installed updates

Turn Windows features on or off

To uninstall a program, select it from the list and then click Uninstall, Change or Repair.

Name	Publisher	Installed On	Size	Version
Microsoft ODBC Driver 11 for SQL Server	Microsoft Corporation	18/04/2020	9.23 MB	12.0.2000.8
Microsoft ODBC Driver 13 for SQL Server	Microsoft Corporation	18/04/2020	8.50 MB	14.0.1000.169
Microsoft Office Enterprise 2007	Microsoft Corporation	06/04/2020	11.1 MB	12.0.4518.1014
Microsoft Report Viewer 2014 Runtime	Microsoft Corporation	18/04/2020	32.5 MB	12.0.2000.8
Microsoft Silverlight	Microsoft Corporation	07/05/2020	3.11 MB	1.0.20926.0
Microsoft SQL Server 2005	Microsoft Corporation	07/05/2020		
Microsoft SQL Server 2005 Backward compatibility	Microsoft Corporation	23/07/2020	45.7 MB	8.05.2309
Microsoft SQL Server 2008 R2	Microsoft Corporation	20/09/2020		
Microsoft SQL Server 2008 R2 (64-bit)	Microsoft Corporation	17/10/2019		
Microsoft SQL Server 2008 R2 Management Objects	Microsoft Corporation	18/04/2020	1.61 MB	10.51.2500.0
Microsoft SQL Server 2008 R2 Native Client	Microsoft Corporation	06/08/2020	18.3 MB	10.52.4000.0
Microsoft SQL Server 2008 R2 Setup (English)	Microsoft Corporation	20/09/2020	26.9 MB	10.52.4000.0
Microsoft SQL Server 2008 Setup Support Files	Microsoft Corporation	18/04/2020	43.3 MB	10.3.5500.0
Microsoft SQL Server 2012 Native Client	Microsoft Corporation	18/04/2020	13.3 MB	11.3.6540.0
Microsoft SQL Server 2014	Microsoft Corporation	02/04/2020		
Microsoft SQL Server 2014 Management Objects	Microsoft Corporation	15/04/2020	14.3 MB	12.0.2000.8
Microsoft SQL Server 2014 Policies	Microsoft Corporation	18/04/2020	2.01 MB	12.0.2000.8
Microsoft SQL Server 2014 Setup (English)	Microsoft Corporation	18/04/2020	77.7 MB	12.0.2000.8
Microsoft SQL Server 2014 Transact-SQL Compiler Ser...	Microsoft Corporation	18/04/2020	281 MB	12.0.2000.8
Microsoft SQL Server 2014 Transact-SQL ScriptDom	Microsoft Corporation	18/04/2020	9.26 MB	12.0.2000.8
Microsoft SQL Server 2017	Microsoft Corporation	18/04/2020		
Microsoft SQL Server 2017 Policies	Microsoft Corporation	18/04/2020	1.03 MB	14.0.1000.169
Microsoft SQL Server 2017 T-SQL Language Service	Microsoft Corporation	18/04/2020	8.49 MB	14.0.17277.0
Microsoft SQL Server Compact 3.5 Design Tools ENU	Microsoft Corporation	07/05/2020	8.40 MB	3.5.5386.0
Microsoft SQL Server Compact 3.5 ENU	Microsoft Corporation	07/05/2020	2.69 MB	3.5.5386.0

Right Click Microsoft SQL Server 2008 R2

And Select Uninstall/Change


SQL Server 2008 R2

Add
Choose this option to add features.


Repair
Choose this option to repair the product.

Remove
Choose this option to remove the product.


Click on Remove Option.




Click OK to Proceed.




Select WELSQL2008R2 from list and Click next




Click on Select all and Click Next




Click Next to Proceed.



Click Remove




Wait for Completion..




Once completed Click Close Button.

Removing Application Files.


In desktop Icons, Right Click on WELTRAUM PMS or WELTRAUM RMS.



See the Location of Main folder (here C:\ WELTRAUM)



Delete the "WELTRAUM" folder



If asking file access denied message Click on “Do this for all current Items” and Click Continue.

Remove all WELTRAUM SOFTWARE icons from desktop.

Uninstallation Process Completed